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Bio-Genesis Technology

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RNA 純化淺談	P.1 - P.3
From 細胞、組織、細菌	
▶ RNA-Bee	P.4
▶ MasterPure™ Complete DNA and RNA Purification Kit	P.5 - P.6
▶ NucleoSpin® RNA II	P.7 - P.8
▶ NucleoSpin® RNA L	P.9
From 單一細胞、LCM sample	
▶ NucleoSpin® RNA XS	P.10
▶ ArrayPure™ Nano-scale RNA Purification Kit	P.11
From 酵母菌、菌絲體	
▶ MasterPure™ Yeast RNA Purification Kit	P.12
Small RNA Application	
▶ NucleoSpin® miRNA	P.13 - P.14
From 植物組織	
▶ NucleoSpin® RNA Plant	P.15 - P.16
▶ MasterPure™ Plant RNA Purification Kit	P.17
From FFPE sample	
▶ WaxFree™ RNA	P.18
▶ NucleoSpin® FFPE RNA	P.19 - P.20
▶ QuickExtract™ FFPE RNA Extraction Kit	P.21
快速萃取	
▶ QuickExtract™ FFPE RNA Extraction Kit	P.21
▶ QuickExtract™ RNA Extraction Kit	P.22
mRNA Isolation	
▶ NucleoSpin® mRNA Mini、Midi	P.23
▶ mRNA-ONLY™ Prokaryotic/Eukaryotic mRNA Isolation Kit	P.24
RNA Clean-Up	
▶ NucleoSpin® RNA Clean-Up、RNA Clean-UP XS	P.25 - P.26
其他相關產品	
▶ NucleoSpin Triprep(可純化 RNA/DNA/Protein)	P.27
▶ NucleoSpin RNA/Protein	P.28
▶ NucleoSpin RNA/Protein	P.29

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RIN Value	Most Suitable Application
1 ~ 4	PCR assays with short region of amplification
4.1 ~ 6.9	qRT-PCR applications
7~10	Highly demanding gene array assays

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動物細胞、動物組織、細菌 → RNA Bee、NucleoSpin RNA II or MasterPure RNA purification kit

單一細胞或是微量動物組織 → NucleoSpin RNA XS or ArrayPure Nano-Scale RNA purification kit

植物組織 → NucleoSpin Plant or MasterPure Plant purification kit

FFPE sample → NucleoSpin FFPE or QuickExtract FFPE RNA extraction kit

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miRNA assay → NucleoSpin miRNA

快速篩檢 → QuickExtract RNA extraction kit

RNA Clean-Up → NucleoSpin RNA Clean-Up or RNA Clean-Up XS

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
Targets	Samples Extracted
Human & mammalian RNA & genomic DNA	Plasma
HIV	Whole blood
HCV	Guthrie cards
RSV	Buccal cells
HPV	Liver
Yeast	Mouse tail
Enterovirus	Kidney
<i>B. pertussis</i>	Saliva
<i>E. coli</i>	Urine
<i>M. tuberculosis</i>	Sputum
Streptococcus mutans	Tissue-culture cell lines
Soy (tofu)	Cervical cells
Maize	Paraffin-embedded tissues
Insect tissues	

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1, 2: + MDB
3, 4: - MDB
5, 6: product Q
7, 8: + control
9: - control
10: marker


Fig.5 rDNase 搭配 MDB (membrane desalting buffer)能完全移除 genomic DNA

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Product Information

Technology	Silica – membrane technology
Format	Mini spin columns
Sample material	< 5×10^6 cultured cells, < 10^9 bacterial cells < 10^8 yeast cells, < 30 mg tissue
Fragment size	> 200 b
Typical yield	14 μ g from 10^6 HeLa cells 70 μ g from 10^9 bacterial cells
A260/280	1.9-2.1
Typical RIN (RNA integrity number)	> 9
Elution volume	40-100 μ l
Preparation time	30 min/6 preps
Binding capacity	200 μ g

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A: RNA isolated from 100mg (1~4) and 200mg (5~8) pig kidney.

Lane 1,2,5,6: NucleoSpin RNA L


Lane 3,4,7,8: Competitor Q

B: RNA isolated from 200mg pig liver sample.

Lane 1,2,3: NucleoSpin RNA L

Lane 4: Blank

Lane 5,6,7: Competitor Q

Product Information

Technology	Silica – membrane technology
Format	Midi spin columns
Sample material	< 5 x 10 ⁷ cultured cells, < 200 mg tissue < 10 ¹⁰ bacterial cells, up to 3 x 10 ⁸ yeast cells
Fragment size	200b ~ 20Kb
Typical yield	180 µg from 10 ⁷ HeLa cells 620 µg from 4 x 10 ⁷ HeLa cells
A260/280	1.9-2.1
Typical RIN (RNA integrity number)	> 9
Elution volume	500 µl
Preparation time	80 min/4 preps
Binding capacity	700 µg

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

Tissue or Cell Source	Amount of Starting Material	NucleoSpin RNA XS Sample No. 1 (ng/μl)	NucleoSpin RNA XS Sample No. 2 (ng/μl)	Competitor Q Sample No. 1 (ng/μl)	Competitor Q Sample No. 2 (ng/μl)
Mouse Kidney	0.5 mg	157	108	96	63
Mouse Kidney	0.1 mg	27	25	17	18
HeLa Cells	1,000 cells	1.3	1.7	0.08	0.06

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Product Information

Technology	Silica – membrane technology
Format	XS spin columns
Sample material	Tissue <5mg 1-10 ⁵ cultured cells
Fragment size	>200b
Typical yield	10 ² HeLa cells: 0.1~1.5ng ; 10 ³ HeLa cells: 10~15ng 10 ⁴ HeLa cells: 100~150ng ; 10 ⁵ HeLa cells: 1000~1500ng
A260/280	1.9-2.1
Typical RIN (RNA integrity number)	> 9 (Depending on sample quality)
Elution volume	5-20 μl
Preparation time	40 min/12 preps
Binding capacity	90 μg

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
產品內容: Nano-scale Lysis Solution, 2X Nano-scale Lysis Solution, MPC Protein Precipitation Reagent, Proteinase K, RNase-Free DNase I, RiboGuard™ RNase Inhibitor, 1X DNase Buffer, TE Buffer.

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
Quality of obtained using the MasterPure™ Yeast RNA Purification Kit.

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Microarray analysis of yeast RNA purified with the MasterPure™ Yeast RNA Purification Kit.

右圖是 *S.cerevisiae* 經 MasterPure™ Yeast RNA Purification Kit 萃取 RNA 並經過 reverse transcription、labeling 並 hybridized to *S. cerevisiae* DNA microarray。

Ordering Information

產品名稱	包裝	產品編號
MasterPure™ Yeast RNA Purification Kit	10 RNA purifications	MPY83010
	100 RNA purifications	MPY83100

產品內容: Extraction Reagent for RNA, MPC Protein Precipitation Reagent, TE Buffer (in 100-reaction kit only), Proteinase K, RNase-Free DNase I, RiboGuard™ RNase Inhibitor, 10X DNase Buffer, 2X T & C Lysis Solution.

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Product Information

Technology	Silica – membrane technology
Format	Mini spin columns
Sample material	< 10 ⁷ culture cells, < 30mg human/animal tissue < 50mg plant tissue, < 150ul reaction mixture
Fragment size for Small RNA	18 - 200 bases
Fragment size for Large RNA	>200 bases
Typical yield	10ug small RNA, 95ug large RNA from 10 ⁷ HeLa cells
Binding capacity	200µg
Elution volume	30 - 100 µl
Preparation time	< 45 min (6 preps human/animal tissue, small and large RNA) < 35 min (6 preps human/animal tissue, small RNA)

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Plant RNA Extraction


Fig.1 NucleoSpin RNA Plant 能移除所有 RT-PCR inhibitor

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Lane 5: rye,

M: Marker.

Lanes	Average yield
1 – 3: wheat MN	24 µg
4 – 6: wheat competitor Q	18 µg
7 – 9: corn MN	36 µg
10 – 12: corn competitor Q	32 µg

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
Lanes:	Average Yield
1-3: poppy seed lysed with buffer RA1	0.84ug
4-6: poppy seed lysed with buffer RAP	7ug

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Product Information

Technology	Silica – membrane technology
Format	Mini spin columns
Sample material	1-100mg tissue
Fragment size	200 b - 20kb
Typical yield	3-70ug from 100mg plant material
A260/280	1.9-2.1
Elution volume	60 µl
Preparation time	30 min/6 preps
Binding capacity	200 µg

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Source of RNA	Yield (µg)
Alfalfa sprouts	36.1
Banana leaves	6.2
Grape leaves	9 - 25.9
Maize leaves	17.7 - 55.7
Maize seedlings	28 - 37 (leaf) 35 (stem) 19 (root)
Maize roots (mature)	18.6
Pine needles	5.2
Raspberry leaves	15.7
Strawberry leaves	12.4
Strawberry roots	46

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
-RNA Extraction from paraffin sample

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
50 reactions

WR-50


WaxFree RNA extraction from lung tissue FFPE section

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
- Lane 1: 100bp ladder;
- Lane 2: NTC
- Lane 3: gDNA control;
- Lane 4: WaxFree RNA extract
- Lane 5: WaxFree RNA extract without tissue
- Lane 6: Commercial total RNA (positive control)

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
	MN	Q	R
Working steps:	36	43	54
Centrifugation steps:	9	9	15
Preparation time:	70 min	75 min	251
Ct-values**:	21.5	26.0	24.0

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Product Information

Technology	Silica – membrane technology
Format	Mini spin columns - XS design
Sample material	Up to 7 sections (10um)of 250mm ² or 15mg paraffin
Typical yield	Depend on amount and quality of sample
A260/280	1.9 - 2.1
Typical RIN	2 - 6 (Strongly depend on sample quality)
Elution volume	5 - 30 µl
Preparation time	~70 min/6 preps
Binding capacity	60 µg

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Lane M, 100-bp DNA ladder;
Lane 1, a 116-bp region of RYR1;
Lane 2, a 162-bp region of ACTA;
Lane 3, a 166-bp region of RSP18;
Lane 4, a 166-bp region of ENO3;
Lane 5, a 226-bp region of GAPDH;
Lane 6, a 232-bp region of OAZ1;
Lane 7, a 307-bp region of ACTB;
Lane 8, a 308-bp region of TNF.

廠牌	操作時間	Tissue Size	Toxic Chemicals	Special requirement	Tube change
Epicentre FFPE DNA Extraction Kit	32 min ★	1 section ★	No ★	No ★	No ★
Q牌 FFPE DNA Extraction Kit	5 hr to Overnight	3 section	Xylene	Column and EtOH ppt.	Yes
S牌 FFPE DNA Extraction Kit	Overnight	5 section	No	Filter Column and EtOH ppt.	Yes
P牌 FFPE DNA Extraction Kit	Overnight	1 mg tissue	No	Resin and DNA IQ System	Yes
G牌 FFPE DNA Extraction Kit	7 hr to overnight	3 section	xylene	EtOH ppt.	yes

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產品內容：QuickExtract™ FFPE RNA Extraction Solution, DNase buffer 1, RNase-Free DNase I, Stop Solution.

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- Lane M, 100 bp ladder;
- Lanes 1 and 2, QuickExtract™ RNA Extraction Kit;
- Lanes 3 and 4, kit from Competitor 1;
- Lanes 5 and 6, kit from Competitor 2.

Ordering Information

產品名稱	包裝	產品編號
QuickExtract™ RNA Extraction kit	5ml (50 reactions)	QER09015
	50ml (500 reactions)	QER090150


產品內容：QuickExtract™ RNA Extraction Solution, DNase buffer 1, RNase-Free DNase I, Stop Solution, RiboGuard™ RNase Inhibitor.

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
產品內容: Terminator™ 5'-Phosphate-Dependent Exonuclease, mRNA-ONLY™ Prokaryotic 10X Reaction Buffer, RiboGuard™ RNase Inhibitor, mRNA-ONLY™ Stop Solution, 5 M LiCl Solution, RNase-Free Water.

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
Product Information

	NucleoSpin® RNA Clean-up	NucleoSpin® RNA Clean-up XS
Technology	Silica – membrane technology	Silica – membrane technology
Format	Mini spin columns	XS spin columns
Sample material	Up to 200µl	Up to 300µl
Fragment Size	>200b	>200b
Typical Recover	85 ~ 95%	85 ~ 95%
A260/280	1.9 - 2.1	1.9 - 2.1
Elution Volume	40 - 100 µl	5 - 20µl
Preparation Time	20min / 6 preps	20min / 6 preps
Binding Capacity	200 µg	90 µg

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Technology	Silica-membrane technology
Format	Mini spin columns
Sample material	≤ 5 x 10 ⁶ cells* ≤ 30 mg human/animal tissue* ≤ 100 mg plant tissue*
DNA yield	≤ 6 µg
RNA yield	≤ 70 µg
Protein yield	≤ 1200 µg
A ₂₆₀ /A ₂₈₀ DNA	1.7-1.9
A ₂₆₀ /A ₂₈₀ RNA	1.9-2.1
Typical RIN (RNA integrity number)	Depending on sample quality: good quality samples (e.g. fresh HeLa cells or mouse liver) typically result in RIN >9
Elution volume DNA	100 µl
Elution volume RNA	40-120 µl
Resolubilization volume protein	10-100 µl
Preparation time RNA and DNA	45 min/6 preps
Preparation time protein	35 min/6 preps
Binding capacity (RNA)**	200 µg

Sample Material	Sample Amount
Human cells (HeLa)	10 ⁶ cells
Mouse liver	3 mg
Fish (Zebrafish)	30 mg
Plant root (garden cress)	100 mg
Plant leave (garden cress)	100 mg
Yeast	10 ⁸ cells
Bacteria	10 ⁹ cells

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1. Lamda DNA marker HindIII
2. HeLa
3. Mouse Liver
4. Zebrafish
5. garden cress root
6. garden cress leave
7. Yeast
8. Bacteria

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1. Lamda DNA marker HindIII
2. HeLa
3. Mouse Liver
4. Zebrafish
5. garden cress root
6. garden cress leave
7. Bacteria
8. Yeast

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NucleoSpin RNA/Protein 簡易操作流程




Fig1. Inducible Knock-down of Hint1 in MCF-7 Tet-On cells

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

⊕: positive control ⊖: negative control

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Bio-Genesis Technology Inc.