PCR & RT-PCR 專刊

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FailSafe™ PCR System

MasterAmp™ Extra-Long PCR Kit

TAQXpedite PCR System

Q-Amp™ 2x HotStart PCR Master Mix

Q-Amp™ 2x ScreeningFire Taq Master Mix

2x pfu PCR Master Mix

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

1.25 ng human gDNA or cDNA was amplified using a gene-specific primer and Q-Amp HotStart PCR Master mix.

The amplicons obtained were:

Lane 1 : 285bp; obtained from human cDNA and amplified using actin primer.

Lane 2 : 564bp; obtained from human cDNA and amplified using actin primer.

Lane 3: 1020bp; obtained from human gDNA and amplified using Glu primer.

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
Product information:

#DCB-001, 1mL

2x ScreeningFire Taq Master Mix testing with human genomic DNA

Template: Human gDNA from blood

Primer set: Glutathione Peroxidase, 1020bp

M: Bio-100 Ladder

Lane 1: 20ng; Lane 2: 10ng

Lane 3: 1ng; Lane 4: 100pg

Lane 5: 10pg; Lane 6: 1pg

Fast protocol

2% Agarose gel in 0.5X TBE

Buffer

M: 100bp marker

Lane1: TOBACCO genomic DNA (600bp products)

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

Lane 1: No-RT control

Lane 2: PCR products from cDNA synthesized by MMLV HP RT


Lane M: 100bp ladder

Products Information

產品名稱	Cat#	Size
MMLV High Performance Reverse Transcriptase 內含 MMLV HP RTase、10x reaction buffer、DTT	RT80110K	10,000 Units (200 U/μl)
	RT80125K	25,000 Units (200 U/μl)
MMLV Reverse Transcriptase 1st-Strand cDNA Synthesis Kit 內含 MMLV RTase、10x reaction buffer、RNase inhibitor、dNTP pre-mix、DTT、Oligo(dT)21 primer、Random 9-mer primer、RNase-Free water	MM070110	10 Reactions
	MM070150	50 Reactions

Q-Amp MMuLV cDNA Synthesis Kit (random primer or oligo dT)

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1. 內含 Random nonainer primer 或 Oligo d(T)
2. 使用 MMuLV reverse transcriptase
3. 只要 1pg total RNA 即可反應 (Tested 1pg of total RNA)

Amplification of the Zebrafish Positive Control Template from the Q-Amp MMuLV cDNA synthesis by two-step RT-PCR. Quantities of 1μg RNA template were used for 20μl reverse transcription reactions. 100ng to 0.01ng of the cDNA SYBR Green qPCR Master Mix and bb3 gene-specific primers were used in a 10μl qPCR reaction with illumina Eco qPCR system (A). The qPCR product was also analyzed by 1% agarose gel (B).

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
Lane M: DNA Marker

Lane 1: MasterAmp with one-step protocol

Lane 2: RT-PCR with MMLV & Taq with two-step procedure

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Pathogen name	Catalogue number	Pathogen genome type
Dengue Virus subtypes 1,2 3 and 4	Path-Deng	RNA
Escherichia coli 0157:H7	Path-e.coli-0157	DNA
Escherichia coli_O104:H4	Path-E.coli_0104:H4	DNA
Hepatitis B Virus	Path-HBV	DNA
Hepatitis C Virus	Path-HCV	RNA
Herpes Simplex Virus 1 and 2	Path-HSV1&2	DNA
Human Immunodeficiency Virus type 1	Path-HIV1	RNA
Human Immunodeficiency Virus type 2	Path-HIV2	RNA
Human influenza A virus subtype (H1)	Path-H1	RNA
Human influenza A virus subtype (H3)	Path-H3	RNA
Human influenza B virus	Path-fluB	RNA
Human Coronavirus 2012	Path-HCoV_2012	RNA
Mycobacterium Tuberculosis	Path-TB	DNA
Mycobacterium tuberculosis complex	Path-M.tuberculosis_complex	DNA
Influenza A virus subtype H5N1 (avian influenza)	Path-H5N1	RNA
Swine H1N1 Influenza Human Pandemic Strain	Path-H1N1-swine	RNA
Vibrio all species	Path-Vibrio_spp	DNA
Vibrio cholerae all subtypes	Path-V.cholerae_subsp	DNA
Vibrio cholerae toxogenic subtypes	Path-V.cholerae	DNA

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1. Primer and probe mix (150 reactions)
2. Reverse Transcription, target specific primers (RNA genome viruses only)
3. Copy number standard curve (sufficient for 4 standard curves)
4. RNase/DNase free water

genesig 高級套組內容物 (Advanced kit contents)

1. Primer and probe mix (150 reactions)
2. Reverse Transcription, target specific primers (RNA genome viruses only)
3. Copy number standard curve (sufficient for 4 standard curves)
4. Internal extraction control - Read through VIC channel* (150 tests)
5. Endogenous ACTB control (150 tests)
6. RNase/DNase free water

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
Product Name	Cat#	Detection chemistry	Model		
Precision OneStep™ qRT-PCR Mastermix	OneStep	TaqMan probe	Illumina	Eco	
			Roche	LightCycler 480, LightCycler® Nano	
	OneStep-SY	SYBR Green	Bio-Rad	MJ Chromo4, CFX96, MiniOpticon, CFX384	
			Thermo	PikoReal	
			Qiagen (Corbett)	Roto-Gene™ 3000 & 6000, Roto-Gene™ Q	
			Cepheid	SmartCycler™	
Eppendorf	Mastercycler® ep realplex				
Precision OneStep™ qRT-PCR Mastermix with LOW ROX	OneStep-LR	TaqMan probe	ABI	7500, 7500Fast, ViiA	
	OneStep-LR-SY	SYBR Green			
Precision OneStep™ qRT-PCR Mastermix with ROX	OneStep-R	TaqMan probe	ABI	7000, 7300, 7700, 7900, 7900HT, StepOne, StepOne plus	
	OneStep-R-SY	SYBR Green			
Precision OneStep™ qRT-PCR Mastermix for the Bio-Rad iCycler	OneStep -iC	TaqMan probe	Bio-Rad	iCycler, IQ4 and IQ5 platforms	
	OneStep -iC-SY	SYBR Green			
Precision OneStep™ qRT-PCR Mastermix for Stratagene hardware	OneStep -MX	TaqMan probe	Stratagene	MX platforms	
	OneStep -MX-SY	SYBR Green			
Precision OneStep™ qRT-PCR Mastermix for capillary lightcyclers	OneStep -CL	TaqMan probe	Roche	Capillary Lightcycler 1.0, 1.5, 2.0	
	OneStep -CL-SY	SYBR Green			

Precision FAST Mastermix and Precision Mastermix for hydrolysis probes

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Detection of 18S Ribosomal RNA using Precision™ MasterMix with ROX for the ABI7000

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