## Transfection reagent

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#### Polyplus

jetPRIME (multifunction transfection)	4
jetPEI(DNA transfection)	5
特殊細胞 DNA transfection reagent	6
in vivo-jetPEI (in vivo transfection)	7
INTERFETin (siRNA & miRNA transfection)	9
jetMESSENGER (mRNA transfection)	9
FectoPRO (bioproduction)	9

#### Mirus

TransIT-X2 (DNA & siRNA transfection)	10
TransIT-2020 (DNA transfection)	11
TransIT-LT1 (DNA transfection)	11
特殊細胞 DNA transfection reagent	12
TransIT-Insect (昆蟲細胞)	13
Ingenio-Electroporation Product (電穿孔試劑)	14

**Polyplus**<sup>+</sup>  
transfection®

Delivery by

**Mirus**

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服務信箱: [service@biogenesis.com.tw](mailto:service@biogenesis.com.tw)

## 細胞轉染 (Transfection)

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Cell types	Cell lines	Description	Transfection efficiency
Epithelial	B16-F10	Murine melanoma	70-80%
	BNL-C12	Murine normal embryonic hepatocyte	50-60%
	CaCO2	Human colon carcinoma epithelial	20%
	CHO-K1	Chinese hamster ovary	70%
	HCT-116	Human colon carcinoma	70%
	HeLa	Human cervix epitheloid carcinoma	70-90%
	HepG2	Human hepatocarcinoma	50-70%
	Huh-7	Human hepatocarcinoma	30-50%
	MCF-7	Human breast adenocarcinoma	50%
	MCF-10A	Human breast adenocarcinoma	40-50%
Fibroblast	MDCK	Canine kidney epithelial	20%
	PC-3	Human prostate carcinoma	70%
	Vero	African green monkey kidney	50%
	COS-7	African green monkey kidney	60-80%
Myeloblast	HEK-293	Human embryonic kidney fibroblast	80-90%
	MRC-5	Human lung fibroblast	50%
	NIH-3T3	Murine embryonic fibroblast	50-70%
Myeloblast	Raw 264.7	Murine monocyte/macrophage	40-50%
Myoblast	C2C12	Murine myoblast	70-90%
Neuronal	SH-SY5Y	Human neuroblastoma	70-80%
Primary Hepatocytes		Human primary hepatocyte cell	20-30%
Primary Melanocytes		Human primary melanocyte cell	40-50%

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## Efficient in a wide range of cell types

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### Tail vein injection

**Nucleic acid:** 50 µg  
**in vivo-jetPEI™:** 5-8 µl  
**N/P ratio:** 5-8  
**Injection volume:** 200-400 µl, 5% glucose  
**Method:** The mouse is placed in a restrainer and 70% ethanol is applied on the tail in order to slightly swell the vein. Complexes in solution are injected into the tail vein over 10 sec, using a ½ inch 26 gauge needle and a 1 ml syringe.


Tissue	Relative gene transfer efficiency
Intestine	~1
Pancreas	~10
Ovary	~100
Brain	~10
Spleen	~100
Kidney	~100
Heart	~100
Liver	~100
Lung	~1000

### Intracerebral injection (stereotaxic)

**siRNA:** 0.1 µg/µl using **jetSI™ 10 mM** (see datasheet)  
**Injection volume:** 1-4 µl  
**Method:** Single injection into either lateral ventricle or stereotaxic

**DNA:** 1 µg (for 8-12 week-old mice)  
**in vivo-jetPEI™:** 0.12 µl  
**N/P ratio:** 6  
**Injection volume:** 5 µl, 5% glucose  
**Method:** Single injection (5 µl) into either lateral ventricle (0.2 mm posterior to the bregma line, 1.1 mm lateral, and 2.2 mm deep from the pial surface) to pentobarbital anesthetized mice (65 mg/kg).

Example of transfected cells expressing the β-galactosidase and found in the anterior subventricular zone (1 week after intraventricular injection of pCMVlacZ).  
 Courtesy B. Demeneix

### Intraperitoneal injection

**Nucleic acid:** 100 µg  
**in vivo-jetPEI™:** 12-16 µl  
**N/P ratio:** 6-8  
**Injection volume:** 400 µl to 1 ml, 5% glucose  
**Method:** Complexes in solution are injected into the peritoneal cavity over 10 sec, using a ½ inch 26 gauge needle and a 1 ml syringe.

Tissue	Relative gene transfer efficiency
Diaphragm	~1000
Uterus	~1000
Salivary gland	~10
Intestine	~10
Muscle	~10
Stomach	~1000
Ovary	~1000
Pancreas	~1000
Brain	~1000
Spleen	~100
Kidney	~100
Liver	~100
Lung	~100

### Intratumoral injection

**Nucleic acid:** 10-20 µg  
**in vivo-jetPEI™:** 1.23-2 µl  
**N/P ratio:** 6-8  
**Injection volume:** 50-100 µl, 5% glucose  
**Method:** For implanted subcutaneous tumors (size > 5 mm<sup>3</sup>), perform multiple injections of 10-20 µl complexes at different sites of the tumor to avoid reflux.

Injection)

page 51)


axial injection.

## Retro-orbital injection

**Nucleic acid:** 40 µg  
**in vivo-jetPEI™:** 6.4 µl  
**N/P ratio:** 8  
**Injection volume:** 200-400 µl, 5% glucose

**Method:** The tip of a 27 g hypodermic needle is introduced carefully in front of the eye. Follow the edge of the orbit down until feeling the needle tip at the base beneath the eye. Inject complexes in solution within 2 sec. If performed carefully, there will be little or no bleeding. The capillary nexus will take up the injected solution rapidly.

Tissue distribution of luciferase transgene expression 24 h following retro-orbital injection.

## Nasal instillation for trachea and lung delivery

**Nucleic acid:** 20 µg  
**in vivo-jetPEI™:** 2.4-3.2 µl  
**N/P ratio:** 6-8  
**Injection volume:** 50-100 µl, 5% glucose

**Method:** Mice are held supine at an angle of 45° with pressure applied to the lower mandible to immobilize the tongue and prevent swallowing. Complexes in solution are then introduced to the nasal planum using a micropipet.

## Subcutaneous injection

**Nucleic acid:** 3-5 µg  
**in vivo-jetPEI™:** 0.3-0.7 µl  
**N/P ratio:** 5-7  
**Injection volume:** 10 µl, 5% glucose

**Method:** Mice are restrained and complexes are injected subcutaneously in the region of interest.

Discovery	Preclinical study	Phase I	Phase II	Phase III
BioCancell, Israel Plasmid DNA Bladder cancer gene therapy		2006/07	2008/13	
Genetic Immunity, Hungary Plasmid DNA HIV immune therapy		2005	2008/13	
CHU- Toulouse, Rangueil Hospital, France Plasmid DNA Pancreas cancer gene therapy		2011/13		
Senesco Technologies, USA Plasmid DNA and siRNA Multiple myeloma therapy		2011/13		
Ottawa Hospital Research Institute, Canada Ex vivo, plasmid DNA Acute myocardial infarction gene therapy	2013			
BiOncoTech, Spain Long dsRNA Melanoma immunotherapy	2011/13			
Trinity Coll. Dublin, Ireland siRNA Modulation of blood-brain barrier	2012			
Bonn Univ., Germany siRNA Melanoma therapy	2013			

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### INTERFERin 訂購資訊

Catalog Number	Amount of reagent
409-01	0.1 ml
409-10	1 ml
409-50	5 x 1 ml

## jetMESSENGER


### mRNA transfection reagent for hard to transfect cells

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Catalog Number	Amount of reagent	Amount of buffer
150-01	0.1 ml	10 ml
150-07	0.75 ml	60 ml
150-15	1.5 ml	2 x 60 ml

Cell type	Cell name	% transfection efficiency
Brain cells	U-87 MG	>90%
	Cortex neurons	40-60%
Epithelial cells	Caco-2	>90%
	MCF-7	>90%
Fibroblasts	MDCK	40-60%
	BJ	>90%
	Primary fibroblasts	70-90%
Hepatocytes	IMR-90	>90%
	MEF	60-80%
Lymphocytes	Hep G2	>90%
	Jurkat	40-60%
Macrophages	K562	40-60%
	Macrophages	50-70%
Monocytes	Primary monocytes	10-30%
	THP-1	80-90%
Stem cells	hMSC	>90%
	mES	50-70%

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
# TransIT®-BrCa Transfection Reagent

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
# TransIT®-Jurkat Transfection Reagent

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# TransIT®-TKO Transfection Reagent TransIT®-siQUEST Transfection Reagent

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Cell Line (Source)	Endogenous Transcript	Knockdown Efficiency
BNL CL.2 (mouse liver)	MAPK1	80%
	MAPK3	83%
HeLa (human cervix)	Lamin A/C	80%
	GAPDH	80%
Hepa1c1c7 (mouse liver)	MAPK1	80%
	MAPK3	75%
	MEK1	75%
	PTEN	80%
HepG2 (human liver)	MAPK1	80%
NIH 3T3-L1	MAPK1	70%
	MAPK3	70%
Secondary Human Astrocytes	Lamin A/C	80%
Primary Mouse Hepatocytes	ABC A1	70%
	Lamin A/C	81%

## TransIT-TKO 訂購資訊

產品編號	容量
MIR 2154	0.4mL
MIR 2150	1mL
MIR 2155	5x 1mL
MIR 2156	10x 1mL

## TransIT-siQUEST 訂購資訊

產品編號	容量
MIR 2114	0.4mL
MIR 2110	1mL
MIR 2115	5x 1mL
MIR 2116	10x 1mL

# TransIT<sup>®</sup>-Insect Transfection Reagent

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