


Plasmid purification



創世紀生技

siRNA synthesis

BDM2017012

Plasmid DNA Purification Kit  買大送小	購買 MN (Mini/Midi/Maxi) 或 GenePharma (SET) 即享加價購優惠	
	加價購一  TransIT-X2 買一送一	加價購二  SYBR Green 5ml NT3500 5mLx2 驚喜價



GenePharma

- High quality siRNA synthesis
- HPLC purification
- Transfection

SET

Catalog	Amount	Modification	Pre-designed siRNA	
A10001	1 set	Pre-designed siRNAs	3 candidates (2 OD for each)	1 OD NC, 1 OD NCFAM and 1 OD PC
A10002	1 set		4 candidates (4 OD for each)	
A10003	1 set	Chemically-modified siRNA	3 candidates (2 OD for each)	
A10004	1 set		4 candidates (4 OD for each)	
A10005	1 set	Fluorescently-labeled RNA	3 candidates (2 OD for each)	
A10006	1 set		4 candidates (4 OD for each)	

Positive Control

p.s. 代客設計與 set 品質保證(至少一對可達 70% 基因抑制)

Catalog	Type	siRNA	Catalog	Type	siRNA
A08001	Positive	LaminA/C	A08006	Positive	Vimentin
A08002	Positive	GFP274	A08007	Positive	P53
A08003	Positive	Luciferase GL2	A08008	Positive	GAPDH
A08004	Positive	MAPK1	A08009	Positive	Cyclophilin B
A08005	Positive	Beta-Actin			

Negative Control

Catalog	Type	siRNA	Catalog	Type	siRNA
A06001	Negative	Control	A07001	Negative	FAM-labeled Control

優惠期限: 20170701-20170930

服務專線: 0800-211667

Email: service@biogenesis.com.tw


台北 TEL: 02-8792-9619

台北 FAX: 02-8792-9620

Website: www.bio-genesis.com.tw






挑選須知



product	NucleoBond Xtra Midi/Plus	NucleoBond Xtra Maxi/Plus	NucleoSpinPlasmid Transfection-grade	NucleoBond Xtra Midi/Plus EF	NucleoBond Xtra Maxi/Plus EF
Catalog	740410 /740412	740414 /740416	740490	740420/740422	740424/740426
Technology	Anion-exchange chromatography		Silica-membrane technology and endotoxin-reduction technology	Anion-exchange chromatography	
Endotoxin level	1-10 EU/µg DNA		< 50 EU/µg DNA	<0.05 EU/µg DNA	
Fragment size	< 300kbp		< 15 kbp	< 300kbp	
Typical yield	400 µg	1000 µg	15-30 µg	400 µg	1000 µg
Preparation time	70 min/prep 30 min/prep (Plus Kit)	75 min/prep 35 min/prep (Plus Kit)	14 min/6 preps	85 min/prep 45 min/prep (Plus Kit)	90 min/prep 50 min/prep (Plus Kit)

數據驗證

買大送小



省時 快速 高存活率



優惠期限:20170701-20170930

服務專線：0800-211667

Email：service@biogenesis.com.tw

台北 TEL：02-8792-9619

台北 FAX：02-8792-9620

Website：www.bio-genesi.com.tw